

Instruction Manual

e-pill[®] MedSmart[®] Voice Locked Automatic Pill Dispenser

Item: 995040



Carefully review and save this setup guide

Components




Locked Automatic Pill Dispenser (x1)
SKU 995040



Medication Tray (x1)
SKU 992040



Power Cord (x1)
Not Available Separately



Dosage Template Set (x1)
Not Available Separately



Key (x2)
SKU 995030



AA Battery (x4)
Not Available Separately

e-pill® MedSmart® Voice Locked Automatic Pill Dispenser



Introduction

The e-pill® MedSmart Voice locked automatic pill dispenser is designed to remind and dispense pills according to the patient's dosing schedule, reminding by means of voice alarm and blinking lights. MedSmart Voice will make the correct dose available at the correct time, keeping all other pills secure via a locked lid.

We recommend that you always keep e-pill MedSmart Voice locked. Two keys are included - one key for the caregiver and one key to be kept as a backup.

Note: e-pill MedSmart Voice is a locked pill dispenser, but it is not rated as being tamperproof. It should NOT be used to administer pain medications or Schedule II/III drugs.

How it Works

At dose time, the alarm will sound and the internal tray will rotate, making the correct dose available. Turn the dispenser upside down to dispense - this will release the correct dose and silence the alarm.

SETUP – STEP 1

Powering On and Opening the Lid

1. Turn the pill dispenser upside down
2. Slide the battery door latch to unlock
3. Remove the battery door
4. Insert the included 4 AA batteries making sure that they are placed in the correct direction
5. e-pill MedSmart Voice is now on
6. Replace the battery door and slide the door latch to its lock position
7. Keep the dispenser upside down and unlock its large round cover. To unlock, insert the key into the silver lock. Turn the key a $\frac{1}{4}$ turn clockwise (to the right). WARNING: Do not turn the key more than a $\frac{1}{4}$ turn when unlocking or locking.
8. Turn MedSmart Voice right side up
9. Plug the power cord into the wall and into the dispenser
10. On the front of the lid, push the PRESS button down until you hear a click. The door is now unlatched. Lift up to open.

Clock - Setting the Clock

1. Press and Hold button A for a few seconds until “:” stops blinking, (located between the hours and minutes) then release
2. Press button B to advance the hours (note AM/PM)
3. Press button C to advance the minutes
4. Press button A to save. The “:” symbol will now blink
5. You’ll be confirmed with “Pill Dispenser Settings Updated”
6. Optional – Set the clock to use the 24-hour time format: e-pill MedSmart Voice is set to use the 12-hour time format by default (AM/PM). To switch to 24-hour format:
 - a. With the current time displayed, press button C to enter the settings menu
 - b. Press button C again to switch between 12-hour (AM/PM) and 24-hour time formats
 - c. Press button A to save

Alarms - Setting and Deleting Alarms

All alarms MUST be set in chronological order starting at 12 AM (midnight)

1. Press button A. "1" will appear next to the alarm clock symbol
2. To set the alarm time press button B for the hour (note AM/PM) and button C for the minutes
3. Press button A to save and to advance to the next alarm
4. For Alarms 2 thru Alarm 6 continue as above
 - a. Button B for the hour, button C for the minutes, button A to save
5. Exit the Alarm Setting Mode by pressing button A repeatedly until the current time is show (occurs after Alarm 6)
6. You'll be confirmed with "Alarm Setup Successfully"

Alarms - Deleting an alarm

1. Press button A until you see the alarm you wish to delete
2. Press button B until --:-- is displayed (occurs after 11:00 p.m. or 23:00).
3. Exit the Alarm Setting Mode by pressing button A repeatedly until the current time is show (occurs after Alarm 6)
4. You'll be confirmed with "Alarm Setup Successfully"

Alarms - Setting the Alarm Volume

3 alarm volumes: Silent (0), Medium Volume (1), Loud Volume (2)

1. With the current time shown, press button C to enter the settings menu
2. Press button B to switch between the volume levels: 0, 1 and 2
3. Press button A to save

SETUP – STEP 2

Filling the Tray

1. Remove the clear tray that is inside the dispenser. This makes filling the tray easier.
2. e-pill MedSmart Voice comes with 9 different dosage templates. Each template is designed for a specific number of alarms and, very importantly, the time of day you are filling the tray. Place the correct template based on your needs into the tray's template slot.
3. With the tray's e-pill sticker at 6 o'clock, fill the tray counter-clockwise. The first dose goes into the compartment immediately to the right of the tray's e-pill sticker.
4. Fill compartments (up to 28 dosages are available) with pills using the dosage template as a guide.
5. Place the tray into the dispenser. Align the tray so that the compartment with the tray's e-pill sticker is directly below the gray button labeled "B".
6. Close the lid, turn the dispenser over, and lock the dispenser

Locking and Unlocking the Dispenser

Lock: When closing the lid, you will NOT hear a click. Please double check that the lid is completely closed before turning the dispenser over to lock it. To lock, insert the key into the silver lock and turn the key counter-clockwise (turn to the left) a $\frac{1}{4}$ turn.

Unlock: To unlock, turn the dispenser upside down, insert the key into the silver lock and turn the key clockwise (turn to the right) a $\frac{1}{4}$ turn.

Note: When locking or unlocking do not use excessive force or turn the key more than a $\frac{1}{4}$ turn.

Dispensing

1. The device will alarm for up to 60 minutes along with blinking lights
2. Slide open the white dispensing door
3. Dispense by turning e-pill MedSmart Voice upside down
4. Dispensing turns off the alarm. You'll be confirmed with "Thank you" after medications are dispensed.

Alarm Duration

At alarm time a voice reminder "It's time for your medications" will sound for up to 60 minutes. During the first 10 minute the alarm will go off once a minute. For minutes 11-59 the alarm will go off every 3 minutes. The lights will blink every 3 seconds for up to 60 minutes.

Missed Dose

If a dose is not dispensed within 60 minutes after the alarm has started, the voice reminder message "Patient may have missed their medications" will sound every minute for up to 4 minutes. After this, the pills will remain available until the next alarm time when the next dose rotates into position, and thus rotating away the missed dose.


Low Battery

We recommend keeping e-pill MedSmart Voice plugged in. However, when the dispenser starts to run low on battery, the low battery symbol will appear. At this time, replace all 4 AA batteries as soon as possible - using only alkaline batteries. The estimate battery life is 6-12 months when not plugged in.

Note: Only use alkaline/non-rechargeable batteries in e-pill MedSmart Voice. The dispenser does not charge batteries.

Power Saving Mode

To save power, e-pill MedSmart Voice will enter stand-by mode by turning off its display after several minutes of no operation. At alarm time, the dispenser will wake up and sound its alarm.



Medication Reminders

Follow us on:



e-pill, LLC
49 Walnut Street, Bldg.4
Wellesley, MA 02481, USA

United States: 1-800-549-0095
International: +1-781-239-2941
www.epill.com

Disclaimer:

Organizing medications correctly, setting and assigning alarm/dose times and determining what medications go into each slot is the responsibility of you, your caregiver, and/or your health care provider. This e-pill Medication Reminder is only a supplement to any memory system that you currently use to remind yourself to take medications. The user maintains full responsibility to administer medication at the appropriate times and to keep the device secure and locked. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including without limitation, damages or injuries resulting from malfunction of the product. e-pill, LLC's liability shall under no circumstances exceed the amount paid by the user for the dispenser.