


e-pill[®] Medication Reminders

Vibrating Alarm Clock

www.epill.com | 1-800-549-0095

Operating Instructions for Shake Awake

KEY FUNCTIONS:

 Selects Vibrating Alarm


 Disable Alarm

 Selects Audible Alarm

 Stop Alarm + Illuminate Display

 Snooze Alarm

 Lock (disables all "Set" buttons)

 Set Time

 Set Alarm

1. TO SET TIME

Slide SET switch to "🕒". Press HOUR/MINUTE to actual time. If your **shake awake** has an AM/PM display, "PM" indicator will appear when afternoon hours are in effect. After setting slide switch to "🔒" (LOCK).

2. TO SET ALARM

Slide SET switch to "🔔" (AL) appears on display.
Press HOUR/MINUTE continuously until desired alarm time appears.
After setting slide switch to "🔒" (LOCK).

3. TO ACTIVATE ALARM

Choose VIBRATING "🔊" or AUDIBLE "🔊" alarm. "🔊" will appear on display.
When alarm activates slide ALARM switch to "🔕" to disable alarm or press "🔔" to stop alarm.
("🔊" will still appear and alarm will repeat every 24 hours.)

4. TO ACTIVATE ALARM WITH SNOOZE FUNCTION (Zzz)

Slide SET switch to "🔔".
Press HOUR/MINUTE continuously until desired alarm time appears.
Slide SET switch to "Zzz" position. (Z Z) will appear on display.
Choose VIBRATING "🔊" or AUDIBLE alarm "🔊".
When SNOOZE alarm activates press "🔔" to stop.
" Z Z " will flash continuously and alarm will repeat every 8 minutes.
To disable SNOOZE alarm return SET switch to "🔒" position. (Z Z) will disappear from display.

(Note: To have SNOOZE alarm repeat every 8 minutes you must disable alarm within 1 minute from start.
Flashing (Z Z) on the display indicates SNOOZE alarm is still activated.)

REMINDER:

shake awake requires 1 AAA and 2 x AA batteries.

Replace AA batteries when vibrating motor becomes weak.

When display becomes dim replace AAA battery.