



e-pill[®] 4 Alarm Vibrating Pocket Pillbox

The e-pill 4 Alarm Vibrating Pocket Pillbox can be set to remind up to 4 times a day by vibrating, beeping or both.

POWERING ON - If the buttons don't beep when pressed or if the display is not on.

1. Remove the back cover of the reminder by sliding the OPEN switch. The switch is located on the bottom edge of the reminder.
2. Remove both battery insulator strips.
 - a. One insulator strip for the AAA battery and one for the round LR44 batteries.
3. Snap the back cover back on.

CLOCK - Set the current time.

1. Slide the LEFT switch to T-SET.
2. Tap the HOUR button to set the correct hour.
 - a. Note: AM does not get displayed / PM does get displayed.
 - b. Optional: Tap the 12H/24H button to switch between 12-hour & 24-hour format.
3. Tap the MINUTE button to set the correct minute.
4. Slide the LEFT switch all the way down to TIME.

ALARMS - Alarms must be set in order starting with the earliest. They will repeat every day.

1. Slide the LEFT switch to the alarm you want to set (1 for the earliest, 2 for the next etc.).
2. Tap the HOUR button to set the alarm's hour.
 - a. Note: AM does not get displayed / PM does get displayed.
3. Tap the MINUTE button to set the alarm's minute.
4. Tap the center round ALARM ON/OFF button to turn on the alarm.
 - a. A pill icon (☐●) above the alarm number indicates the alarm is turned on.
 - b. Important: Each alarm must be turned on individually.
5. Add the medications intended for the set alarm into that alarm's pillbox compartment.
6. Repeat steps 1-5 for each alarm you want to set.
7. Slide the LEFT switch all the way down to TIME after setting and turning on all your alarms.

NOTIFICATION TYPE - The chosen notification type applies to all alarms.

1. Slide the RIGHT switch to your preferred notification type. VIB, BEEP or VIB+BEEP.

AT ALARM TIME - What happens and what to do when the alarm goes off.

1. The alarm will go off for 30 seconds. It will repeat twice if not turned off, after 4 minutes and after 8 minutes.
2. A pill icon (☐●) and the alarm number (1, 2, 3 or 4) set for that time will appear.
 - a. The alarm number matches the numbered pillbox compartment holding the medications for that alarm.
3. Take your medications.
4. Tap the center round ALARM ON/OFF button to silence the alarm.
 - a. The alarm will go off again tomorrow at the same time (no reset required).

BACKLIGHT - Turning on the backlight.

1. Tap any button to turn on the backlight.

BATTERY - How to replace the batteries.

This e-pill uses 4 replaceable batteries. 1 AAA battery powers the alarms and 3 LR44 batteries powers the display. Replace all batteries when the display fades or when the alarm is getting weak. Battery life is roughly 1 year and depends on the numbers of alarms and notification type set.

1. Remove the back cover of the reminder by sliding the OPEN switch. The switch is located on the bottom edge of the reminder.
2. Remove all 4 batteries.
3. Insert 1 new AAA battery.
4. Insert 3 new LR44 batteries. The “+” side (the side with text) will face up.
 - a. This battery type can be found in any pharmacy or grocery store. Look for “LR44” anywhere on the package.
5. Snap the back cover back on.

RESET - How to reset the 4 Alarm Pocket.

1. Use a thin paper clip.
2. Insert the paper clip into the small hole above the word “RESET.” This is located on the back of the reminder between compartments 2 and 4.
3. Tap down on the paperclip. The reminder will vibrate, beep and 12:00 will be displayed.
4. The reminder is now reset.
5. Set your current time and your alarms.



Disclaimer: This e-pill Medication Reminder is only a supplement to any memory system that you currently use to remind yourself to take medications. The user maintains full responsibility to administer medication at the appropriate times. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. e-pill, LLC's liability shall in no circumstances exceed the amounts paid by the user for the reminder.



e-pill® Medication Reminders

e-pill, LLC. 49 Walnut Street Bldg. 4, Wellesley, MA 02481, USA

Phone USA: 1-800-549-0095

Phone International: +1-781-239-2941

www.epill.com