



Medication Reminders

Operating Instructions
e-pill® Multi-Alarm POCKET

SKU: 759030 **UPC:** 837066000701

Thank you for purchasing an
e-pill Medication Reminder

We suggest you to read the following instructions carefully and keep them for future reference.

This easy-to-use e-pill Multi Alarm POCKET reminder can remind you to take your medications up to 6 times per day by a beeping alarm. Alarms can be set on the hour around the clock. Once one or several (up to six) alarms have been set they will repeat every day, no reset required.

Please activate the timer by removing the plastic cover from the battery located in the battery compartment.

www.epill.com

Version 171122

Setting up Your e-pill® Multi-Alarm POCKET

Setting Time



1. Open the clear lid on the timer.
2. Press the "H" Button to set the hour.
3. Press the "M" Button to set the minute.
"H"=Hour
"M" =Minute
4. Press and Hold "Stop" and then "H" to switch between 24 Hour mode (military time) and 12 Hour mode (AM/PM Time).
A as in AM, P as in PM is shown on the display.
Numbers will stop flashing in 5 seconds after you finished setting the time.
5. Close the clear lid.

Setting Daily Alarms 1-6

1. Open the clear lid on the timer.
Top row is AM (morning) alarms (Red) from 1AM to Noon
Bottom row is PM (evening) alarms (Black) from 1PM to Midnight
2. Push down the hour switch on the top row (AM) or bottom row (PM) to turn on each alarm.
3. Alarm bell symbol and a number will appear on the screen for each alarm. You can set up to 6 daily alarms.
4. Close the clear lid.

TIP: The eraser end of a pencil can be used to set the switches.



Red=AM

Black=PM

Example: Two Alarms Set for
8AM and 4PM

Alarms automatically repeat each day at the same time, no reset required.

Stop the Alarm

1. When you hear the alarm sound, take your medications.
2. Then Press the "STOP" Button to stop the Alarm.

You can set up to Six (6) Daily Auto-Repeating Alarms.

Medical Snooze

Alarm will repeat every 5 minutes twice if the alarm is not stopped. Alarm Duration: 15 Seconds

Delete an Alarm

Push up the hour switch to the "Off" Position for the top row (1AM to Noon) or bottom row (1PM to Midnight). The alarm symbol will disappear from the display.

Disclaimer

This e-pill Medication Reminder is only a supplement to any memory system that you are currently using to remind yourself to take medications. The user maintains full responsibility to administer medications at the appropriate times. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. e-pill, LLC's liability shall in no circumstances exceed the amount paid by the user for the timer.

Battery Replacement

This e-pill uses a replaceable battery. Battery life is approximately 2-3 months depending on the number of alarms you set. Battery door is on the right side of the reminder. Simply insert a fresh CR2032 3V Button Cell with the "+" side up and slide it in. Your e-pill Multi-Alarm POCKET should start right up.

e-pill[®] Medication Reminders

e-pill, LLC
49 Walnut Street, Bldg. 4
Wellesley, MA 02481, USA
1-(800)-549-0095
International: +1 (781)-239-2941
Fax +1 (781)-235-3252



www.epill.com