

e-pill[®] Parkinson's 7 Day x 5 Doses per day weekly organizer with countdown interval timer. Black canvas bag, and black timer with lanyard (960620).

Operating Instructions

If you have Parkinson's, managing your daily medication regimen can be a daunting task. As your Parkinson's disease progresses, it's often harder to manage and keep track of your medication schedule and intake. Understanding your medication and keeping to a schedule will help avoid unpleasant "off" periods due to missed doses.

The e-pill Parkinson's weekly pillbox organizer and medication reminder interval timer system may be an ideal solution to help you, "Never forget to take your medications again." This system comes complete with a medication reminder timer and seven pillboxes each with five (5) doses / medication compartments per day in a discreet canvas bag. Organize all your daily medications once a week.

The medication reminder interval timer has a large digital display. The timer has a 12/24-hour clock, interval timer can be set up to 10hr. It's very simple to set and start the timer. Please note you will need to reset the timer by pressing the Start/Stop button.

Never Forget to take your Medications Again!

SETTING your e-pill Interval Timer:

CLOCK – Set the Current Time

- Press 12hr = 24hr button (on the back of the unit) to switch between 12 hour and 24 hour setting (You will notice the AM sign to the left of the clock flash on & off.).
- Press and hold CLOCK button until display disappears (about 3 seconds).
- Press START/STOP to switch between AM or PM.
- Press the number keys "0" through "9" to set the time. Press CLEAR if an incorrect number is selected. Note: for the hours between 1 and 9. '0' must be entered before selecting the hour. Once the hour, minutes, and seconds are entered (**seconds MUST be entered, for example: setting your clock to 8AM , Your display should look like this ^{AM} 8:00 00**). Press CLOCK once to store setting.

COUNT UP

- Press COUNT-UP button and the words "count up" will appear in the right hand corner of the display.



- Press START/STOP to start count-up feature.
- Press START/STOP again to pause counting.
- Press START/STOP again to resume counting.
- Press START/STOP then CLEAR to reset timer to 0:00:00.

TIMER:

- Press TIMER button and the word “Timer” will appear in the display.
- Press the number keys “0” through “9” to desired timing period.

NOTE: to lock the keypad simply move the switch on the side of the timer (near the strap) to 'LOCK'. This will disable the keypad. To unlock the keypad, move the switch to 'UNLOCK'. All buttons are disabled when keypad is locked.)

****For example if you want your timer to go off every 3hrs. Press 3 once and 0 four times****

- Press START/STOP to start countdown.
- Press START/STOP again to pause counting.
- Press START/STOP a third time to resume counting.
- Press START/STOP to stop the alarm sound. The timer will remember the time entered.
- Press CLEAR to reset the timer to 0:00:00.

MEDICATION TIME

BEEPING ALARM – What to do when the ALARM goes OFF

- TAKE your MEDICATIONS
- Press START/STOP to silence alarm after you have taken your medications. Press START/STOP to resume timer.

The alarm will sound for 60-second.

BATTERY – How to replace the Battery

- Your e-pill Timer runs on one (1) replaceable “LR44” battery. The normal battery life is about 12 months.
- The battery door is located on the back of the unit.
- Slide open battery door in direction of arrow.
- Insert one “LR44” battery into compartment. The “+” sign should be facing out.
- Replace the cover by sliding back into place.



Disclaimer: This e-pill Medication Reminder is only a supplement to any memory system that you currently use to remind yourself medications. The user maintains full responsibility to administer medication at the appropriate times. e-pill, LLC assumes no responsibility incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. e-pill, LLC's liability shall in no circumstances exceed the amounts paid by the user for the timer.

