



e-pill® Multi-Alarm POCKET Quick Guide

Medication Reminders

Operating Instructions e-pill® Multi-Alarm POCKET



Please read the instructions carefully and keep them for future reference

Features

- Up to 6 alarms per day
- Alarms can only be set on the hour
- Alarms will repeat every day

Disclaimer

This e-pill Medication Reminder is only a supplement to any memory system that you are currently using to remind yourself to take medications. The user maintains full responsibility to administer medications at the appropriate times. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from

e-pill, LLC
 49 Walnut Street, Bldg. 4
 Wellesley, MA 02481, USA
 1-(800)-549-0095
 International: +1 (781)-239-2941
 Fax: 781-235-3252

Activate


1. Insert 2x AAA
2. Display will turn on

Set Clock

1. Open the clear lid
2. Switch 12h/24h:
Press & hold "STOP" and then tap "H"
A=AM & P=PM
3. Set hour: Press "H"
4. Set minute: Press "M"
5. Time will be saved and will stop flashing 5 seconds after setting the clock

Set Alarms (up to 6)

1. Set alarm: Push down desired switch(s)
 - Red sets alarms: 1AM to Noon
 - Black sets alarms: 1PM to Midnight
2. Alarm bell symbol and a number will appear on the screen for each alarm



Example
 Two alarm set:
 9AM and 3PM

Check Alarms

1. After all alarms are set, tap the "STOP" button to cycle through your set alarms. A=AM & P=PM
2. Close the clear lid

Stop the 15 Second Alarm

1. When you hear the alarm sound, take your medications
2. Then Press the "STOP" Button to stop the Alarm

Medical Snooze

1. If the alarm is not stopped, within 5 minutes, the 15 second alarm will repeat one time
2. During the entirety of the 5-minute medical snooze the corresponding alarm bell will blink

Battery

This reminder uses 2 AAA batteries
 Battery life is approximately 1 year
 The battery door is located on the back of the reminder