



e-pill® Multi-Alarm POCKET Quick Guide

Medication Reminders

Operating Instructions e-pill® Multi-Alarm POCKET

Please read the instructions carefully
and keep them for future reference



Features

Set up to 6 on-the-hour alarms
All alarms will repeat every day

Disclaimer

This e-pill Medication Reminder is only a supplement to any memory system that you are currently using to remind yourself to take medications. The user maintains full responsibility to administer medications at the appropriate times. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. e-pill, LLC's liability shall in no circumstances exceed the amount paid by the user for the timer.

e-pill, LLC
49 Walnut Street, Bldg. 4
Wellesley, MA 02481, USA
1-(800)-549-0095
International: +1 (781)-239-2941
Fax: 781-235-3252

Activate

1. Insert the included 2 AAA batteries
2. The display will turn on

Set Clock

1. Open the clear lid
2. Switch between 12h and 24h mode:
Press & hold "STOP" and then tap "H"
A=am & P=pm
3. Set hour: Press "H"
4. Set minute: Press "M"
5. Time will be saved and will stop flashing
5 seconds after letting go of all buttons

Set Alarms (up to 6)

1. Set alarm: Push down desired switch(s)
 - Red sets am alarms: Midnight to 11am
 - Black sets pm alarms: Noon to 11pm
2. Alarm bell symbol and a number will appear on the screen for each alarm



Example
Two alarm set:
9am and 3pm

Check Alarms

1. Tap the "STOP" button to cycle through set alarms. A=am & P=pm
2. Close the clear lid

Stop the 15 Second Alarm

1. When you hear the alarm sound, take your medications
2. Then press the "STOP" Button to stop the Alarm

Reminds Again

1. If the alarm is not stopped within 5 minutes, the 15 second alarm will repeat one time
2. The corresponding alarm bell will also blink for 5 minute or until the stop button is pressed

Battery

1. This reminder uses 2 AAA batteries
2. Battery life is approximately 1 year
3. The battery door is located on the back of the reminder