



Medication Reminders

e-pill® 6 Alarm Pocket

Please read the instructions carefully and keep them for future reference



Features

Set up to 6 on-the-hour alarms
All alarms will repeat every day

Disclaimer

This e-pill Medication Reminder is only a supplement to any memory system that you are currently using to remind yourself to take medications. The user maintains full responsibility to administer medications at the appropriate times. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. e-pill, LLC's liability shall in no circumstances exceed the amount paid by the user for the timer.

e-pill, LLC

49 Walnut Street, Bldg. 4
Wellesley, MA 02481, USA
1-800-549-0095

International: +1 781-239-2941
Fax: 1-781-235-3252

e-pill® 6 Alarm Pocket Quick Guide

Activate

1. Insert the included 2 AAA batteries
2. The display will turn on

Set Clock

1. Open the clear lid
2. Switch between 12h and 24h mode:
Press & hold **STOP** and then tap **H**
A=am & P=pm
3. Set hour: Tap **H**
4. Set minute: Tap **M**
5. Time will be saved and will stop flashing
5 seconds after letting go of all buttons

Set Alarms (up to 6)

1. Set alarm: Push down desired switch(s)
 - Red sets am alarms: Midnight to 11am
 - Black sets pm alarms: Noon to 11pm
2. Alarm bell symbol and a number will appear on the screen for each alarm

Check Alarms

1. Tap the **STOP** button to cycle through set alarms. A=am & P=pm
2. Close the clear lid

Stop the 15 Second Alarm

1. When you hear the alarm sound, take your medications
2. Then tap the **STOP** Button to stop the alarm

Reminds Again

1. If the alarm is not stopped within 5 minutes, the 15 second alarm will repeat one time
2. The corresponding alarm bell will also blink for 5 minute or until the stop button is pressed

Battery

1. This reminder uses 2 AAA batteries
2. Battery life is approximately 1 year
3. The battery door is located on the back of the reminder