

e-pill Med-Time XL

Automatic Pill Dispenser

www.epill.com **Quick Set-Up Guide** 1-800-549-0095

CLOCK - Set Correct Time

Push button 1 until "CLOCK" is shown on the display. Wait a few seconds. Set clock by pressing button 2 (-) or 3 (+). Push button 1 to store time setting.

DOSES - Set Number of Doses per Day

Push button 1 until "DOSES" is shown on the display. Wait a few seconds until 01, 02 etc. is shown on the display. Set number of daily doses (medication times) by pressing button 2 (-) or 3 (+) until desired number of doses per day is set. Maximum doses per day = 28. Push button 1 to store # of doses per day.

ALARMS - Set Doses Alarm Times

Push button 1 until "ALARMS" is shown on the display. Wait a few seconds. The time for the first alarm is shown. Set desired alarm time by pressing button 2 (-) or 3 (+). If more than one dose per day, set next alarm time by pressing button 1, wait a few seconds and then set next alarm time. Continue until all alarms are set.

Note: Refer to manual for complete use guidelines, instructions and description of features.



**How are
Medications
Dispensed?**

As easy as 1-2-3

- 1. Alarm + Inside Tray Rotates**
- 2. Patient Picks Up Dispenser**
- 3. Right Meds "Pour" Out**

e-pill Med-Time XL



www.epill.com