

# OPERATING INSTRUCTIONS

## e-pill® HÄGER® VIBRATING 8 ALARM MEDICATION REMINDER WATCH

Watch Made in China. Manufactured by HÄGER.

Exclusive Distributor in  
the USA, Canada and Mexico:  
**e-pill®** Medication Reminders



### FEATURES


The e-pill HÄGER Vibrating Medication Reminder Eight (8) Alarm Watch is an innovative device to help improve patient compliance.

- Up to 8 daily alarms
- Choice of vibration, beep or both
- No alarm reset required
- Medical snooze
- Easy-to-read large display with back light

### HOW TO SET THE TIME/ DATE

Your new e-pill HÄGER watch may be in power saving mode (no display). To start your watch press A+B+C+D at the same time. If your watch does not show the correct time please follow these instructions:

- Press & hold button D (SET) for 3 seconds to enter setting mode. The **Hours** will start flashing.
- Tap A (+) or B (-) until the correct **Hour** is reached. Please note the AM/PM indicator in the top left hand corner.
- Tap D move on to **Minutes**. Tap A (+) or B (-) until the correct **Minute** is reached.
- Tap D to move to **Year**. Tap A (+) or B (-) until the correct **Year** is reached.
- Tap D to move on to **Month**. Tap A (+) or B (-) until the correct **Month** is reached (1=January, 2=February ... 12= December).
- Tap D to move on to **Day** of the Month. Tap A (+) or B (-) until the correct **Day** is reached.
- Tab D to leave setting mode.



### TEST ALARM SIGNALS

While in normal time display, press & hold button B (LIGHT) for 5 seconds and the vibrator and beeper will turn on for 2 seconds.

### BACK LIGHT

Tap button B (LIGHT) for Back Light. The Display will be backlit for 3 seconds.

### AM/PM (12 HOUR) or 24 HOUR MODE




Press and hold A button to toggle between AM/PM (12 Hour) or Military Time (24 Hour) mode.

## HOW TO SET DAILY ALARMS

While in normal time mode, Tap button C (ALARM) to go to the alarm you want to set (AL-1 up to AL-8).

- Press & hold button D (SET) for 3 seconds to enter setting mode. The **Hours** will start flashing.
- Tap A (+) or B (-) until the correct **Hour** is reached. Please note the **AM/PM** indicator above in the top left hand corner.
- Tap D to move on to **Minutes**. Tap A (+) or B (-) until the correct **Minute** is reached.
- Tap D to leave setting mode.

Tap A to select alarm mode for each alarm. The possible settings are:

- Sound Alarm On 
- Vibrating Alarm On 
- Both Sound and Vibration Alarm On 
- Alarm OFF (no symbol shown)

To proceed to your next alarm, tap C (ALARM).

## WHEN AN ALARM GOES OFF

The alarm will vibrate, sound or both for 20 seconds, followed by 40 seconds of snooze. This sequence will repeat two (2) times.

### TO STOP THE SNOOZE, TAP ANY BUTTON

**NO Alarm reset required-** Regardless of what you do, alarms will reset at midnight and automatically repeat every day at the same time.

## POWER SAVING FUNCTION

You may want to temporarily shutdown your e-pill HÄGER in order to save the battery (power saving mode). To do this, press all four buttons A+B+C+D at the same time while in normal time display. While in power saving mode, **the clock will stop**, the display will turn off, but the alarm settings will be remained.


To re-start your e-pill alarm watch, press and hold all four buttons A+B+C+D and then set your watch to the correct the Time and Date. Your alarm settings are still retained.

## BATTERY REPLACEMENT

**Battery:** CR2032 3V Button Cell

We suggest you take the watch to a jeweler to have the battery replaced to minimize the risk of mishandling the watch. If you do need to replace the battery yourself, follow these instructions:

- Unscrew the watch bottom cover and open it.
- Unscrew the battery holder.
- Replace the battery with a new 3 Volt Lithium Button Cell (CR2032). The + side should be facing you.
- Screw back the battery holder.
- Use a metal tool to short circuit the two RESET contacts (see illustration below).
- Check the watch display to make sure it is working. If not repeat these steps.
- Screw back the bottom cover with the words at the bottom cover facing up.
- Perform the alarm signal test to make sure vibrator and beeper works.
- Battery life is estimated to be 4 to 6 months but depends on the number of alarms and how long each alarm is sounding.



## WARRANTY

HÄGER covers manufacturing defects for one year from date of purchase. Battery, straps, glass or damage caused by misuse are not covered. During the warranty period, and upon proof of purchase, your watch will be repaired or replaced at the manufacturer's option.

**Contact:** e-pill, LLC (or your distributor)

**Address:** 49 Walnut Street, Bldg. 4  
Wellesley, MA 02481  
USA

**Web site:** [www.epill.com](http://www.epill.com)

**Phone:** (781) 239-2941 or Toll-Free: 1-800-549-0095

**Fax:** (781) 235-3252

### LIABILITY

This watch is designed to help the wearer to remember his/her medication on time. The manufacturer HÄGER takes no responsibility whatsoever for any inconvenience or non-compliance to medications in the event the watch fails to function.

### WARNING

This Operating Instruction is protected by copyright.  
©Copyright e-pill, LLC, 2010  
e-pill® and HÄGER® are trademarks owned by e-pill, LLC and Monico Properties Ltd.