

Operating Instructions

e-pill® MED EZ Keychain Medication Reminder (421131)

e-pill MED EZ Keychain Medication Reminder can remind you to take your daily medications. This Compact, Simple to SET and Easy to USE interval (count-down) timer was designed to be a reminder for medications taken once-a-day or on a 12, 6, 5 or 4 hour interval. It works well as a pill timer for non-chronic medications such as anti-biotics or other Rx drugs taken on a short term basis.

Please note that on this e-pill the battery can be replaced, but we recommend that you take it to a jeweler or place where watch batteries (example Radio Shack) are available. It uses a CR2016.

Never Forget to take your Medications Again!

The interval timer can be set to alarm you:

SCHEDULE:

- Once-A-Day (=Interval 1)
- Every 12 Hours (= Interval 2)
- Every 6 Hours (= Interval 3)
- Every 4 Hours (= Interval 4)
- Every 5 Hours (= Interval 5)

SETTING your e-pill MED EZ Keychain Timer

- Press the front button to open your MED EZ.
- Hold the STOP button down for 10 seconds and the word MED will appear on the display.
- Tap the SET button until you reach the interval you need (see above)
- Press START
- Close the lid



Your e-pill MED EZ Keychain Medication Reminder is now set to remind you to take your medications on the interval you have selected.

MEDICATION TIME How to use your new e-pill Medication Reminder:

ALARM - What to do when the ALARM goes OFF

Alarm will sound for 20 seconds when it is time for you to take your medications. To stop the alarm, open the lid and press the STOP button

- Take your Medications

After taking your medications, press STOP and the MED EZ will continue counting down to the next medication time.

If you need to re-set your MED EZ simply hold down the STOP button for 10 seconds and follow the SETTING instructions above.

Please turn

If medication times are missed, you should follow the directions from your physician or pharmacist before taking any medication out of sequence.



BATTERY

Your e-pill MED EZ Keychain Medication Reminder uses a CR2016 3V Button Cell.

We do not recommend that you try to replace the battery your self, but take it to a jeweler or a Radio Shack.

How to Change the CR2016 Battery:

The battery may be replaced by carefully removing the four (4) screws in the back. Open the case (similar to a watch case) and remove one of the screws holding the battery.

Pry out the battery and replace it with a fresh CR 2016 Button Cell.

Re-assemble your e-pill MED EZ and follow the SETTING instruction.

Disclaimer: This e-pill Medication Reminder is only a supplement to any memory system that you currently use to remind yourself to take medications. The user maintains full responsibility to administer medication at the appropriate times. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. e-pill, LLC's liability shall in no circumstances exceed the amounts paid by the user for the timer.



e-pill® | Medication Reminders

e-pill, LLC, 70 Walnut Street, Wellesley, MA 02481, USA
Phone: 1-800-549-0095 (781-239-8255 U.S & International)

www.epill.com