

e-pill® 4 Alarm Vibrating Pill Box (960622)

This e-pill Medication Reminder can remind you to take your medications up to four (4) times per day by a vibrating / beeping / vibrating + beeping alarm. Alarms repeat automatically every day, no reset required. The four (4) alarms can be set independently of each other. The pill box has four (4) compartments (four doses per day). The easy to read display features a backlight.

Never Forget to take your Medications Again!

SETTING your new e-pill 4 Alarm Vibrating Pill Box:

CLOCK – Set the Current Time

If the display is NOT on, please remove TWO (both) insulator strips located under the pill tray. One insulator is for the clock and the other insulator for the alarm.

- Slide the left side switch to T-SET.
- Press or hold the HOUR button to set the correct HOUR. Please observe the PM indicator.
- Press or hold the MINUTE button to set the correct MINUTES.
- If necessary press the 12H/24H button after setting the time to switch between 12 hour (AM/PM) format and 24 hour (military time) format.
- Slide the switch back to TIME.



ALARMS - Set Alarms (up to 4 daily alarms)

For each alarm slide the left switch to 1, 2, 3 or 4 (Alarm #).

- Slide the left switch to the alarm you want to set.
- Press or hold the HOUR button to set the correct HOUR. Please observe the PM indicator.
- Press or hold the MINUTE button to set the correct MINUTES.
- Press **ALARM ON/OFF** button to turn alarm **ON**. A “pill” symbol must appear with the alarm # (Note. You **MUST** turn on each individual alarm by **PRESSING** the **ALARM ON/OFF** button).
- Slide the switch back to TIME. Repeat for each alarm you want to set.



MODE - Set type of Alarm (VIB / BEEP / VIB + BEEP)

Select how you want to be notified when the alarm goes off. The setting applies to all alarms.

- Slide the right switch and leave it in position. Select VIB (Vibrating Alarm), BEEP (Beeping Alarm) or VIB+BEEP (both).

MEDICATION TIME. How to use your new e-pill Medication Reminder:

ALARM - What to do when the ALARM goes OFF

Alarm will vibrate/sound (or both) for 30 seconds. The alarm will repeat twice - in four (4) minutes and after 8 minutes. While the alarm is vibrating/beeping you will see a flashing “Pill Symbol” and the (1, 2, 3 or 4) and the Scheduled Alarm Time will be shown on the display, directing you to the pill box with the same # on the back of your e-pill device.

- TAKE your MEDICATIONS.
- Press ROUND BUTTON (ALARM ON/OFF) after you have taken your medications.

The alarm will go off again tomorrow at the same time (no reset required).

If medication times are missed, you should follow the directions from your physician or pharmacist before taking any medication out of sequence.

PILLBOX – Daily Medication Tray

Your e-pill comes complete with a Medication Tray (“Pill Box”) containing four compartments. Pill compartments are labeled **1** thru **4** on the back. The medication compartment snaps open and will store your pills when locked.



BATTERY – How to Replace the Batteries

Your e-pill 4 Alarm Vibrating Pill uses four replaceable batteries. One (1) AAA battery is used for the alarm (sound and vibration). Three (3) Button Cell batteries are used for the clock. Replace batteries when the display fades (clock batteries) or when the alarm is getting weak (alarm battery). Normal battery life is 12 months.

To open the battery compartment, slide the switch labeled OPEN in the direction of the arrow, and SLIGHTLY OPEN the top casing. Be careful, if you use too much force to open you may risk breaking the hinges on the back of the top casing. Then lift the back of the casing to separate it from the lower casing. The batteries are now visible.

Remove the old batteries: Place a new fresh AAA (1) battery in the compartment as indicated by the polarity symbols. Place three (3) fresh Button Cell with the positive side up in the battery compartment. Compatible Button Cell batteries:

AG13	280-08 (Citizen)
GPA76	SR44SW (Maxell, SONY, Toshiba, Seiko,)
LR44	9 (Renata)
SR44	S200, A
303 (Duracell, Varta, Ray-O-Vac, Panasonic)	GS14
S06	707

RESET – How to reset the “computer”

Take a paper clip (not a sharp instrument) and insert it in the small hole above the word RESET (located on back of the e-pill 4 alarm vibrating pill box, above the word “Medication” between pill box 2 and pill box 4 - NOT the hole for the sound maker on the right side).

- Press down until 12:00 is displayed. Set the Time and Alarm(s) again.

BACKLIGHT – Turning on the backlight

To turn on the backlight, simply press any button.

Disclaimer: This e-pill Medication Reminder is only a supplement to any memory system that you currently use to remind yourself to take medications. The user maintains full responsibility to administer medication at the appropriate times. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. e-pill, LLC's liability shall in no circumstances exceed the amounts paid by the user for the timer.



e-pill[®] Medication Reminders

e-pill, LLC, 49 Walnut Street Bldg. 4, Wellesley, MA 02481, USA
Phone: 1-800-549-0095 (+1-781-239-2941 U.S & International)

www.epill.com