


This easy to use e-pill Multi-Alarm TimeCap allows you to know when you last took your medication, even without setting any alarms. The automatic LAST OPENED indicator always shows the time and day of week when the TimeCap device was last put on the bottle. Additionally you can set up to 24 daily alarms, on the hour. The e-pill Multi-Alarm TimeCap has been shown to improve patient compliance/ medication adherence in a clinical study.




## SET TIME & DAY OF WEEK

In order to use the features LAST OPENED indicator and setting alarms, you need to set the time and the last day of week.




### 1. Remove the cap from the bottle.

2. Press the MODE button once and the display will start to blink (a).
3. Use the DOWN or UP buttons to change the time. Press and hold either button to move faster (b).
4. When desired time is reached, press the MODE button once and the weekday indicator will start to blink (c).



5. Use the DOWN or UP buttons to change the day of the week (d).
6. When desired weekday is reached, either press the MODE button once or wait 5 seconds to return to the time display (e).




## SET ALARMS



### 1. Remove the cap from the bottle.

2. Press and hold the DOWN and UP buttons at the same time until the display blinks 12 PM (a).
3. Use the DOWN or UP buttons to reach the hour you wish an alarm to sound (b).
4. To set the alarm, press the CHECK button and an alarm symbol will appear in the upper left corner (c).



5. Repeat steps 3 & 4 to set more alarms.
6. To remove an alarm, go to the selected alarm hour and press the CHECK button, and the alarm symbol will disappear.
7. Press and hold the DOWN and UP buttons at the same time, or wait 5 seconds, to store the settings and return to the time display (d).



## HOW TO USE YOUR E-PILL MULTI-ALARM TIMECAP MEDICATION REMINDER


The display on your e-pill Multi-Alarm TimeCap will always record the time that the cap was last opened to remind you of when you last took your medication, even if no alarms are set.



(a)



(b)



(c)

1. **If alarms have been set and an alarm goes off**, your e-pill Multi-Alarm TimeCap will start to beep and the display will start to blink (a).
2. Remove the cap and take your medication.
3. **DO NOT FORGET TO PUT THE CAP BACK ON.**
4. The time and day of week of when you last put the cap back on will automatically be recorded and displayed on the cap (b).

**If an alarm is over due**, the display will stay blinking, showing the last time the cap was opened, to remind you of when you last took your medication. The display will keep blinking until the bottle is opened (c).

## CHECK CURRENT TIME

**NOTE:** The display by default **does not show the current time**. To check the current time, follow these steps.




(a)



(b)



(c)


1. Press and hold the CHECK button, and the display will show the current time and day (a → b). 

2. Let go of the CHECK button for the display to return to the last time you opened the bottle (c)

## CHECK ALARM SETTINGS




(a)




(b)



(c)

1. **DO NOT remove the cap from the bottle.**
2. Press and hold the DOWN and UP buttons at the same time until the display shows 12 PM (a). 

3. Use the DOWN or UP buttons to advance by one hour. The alarm symbol will appear on the display at each hour an alarm is set (b). 
4. Wait 5 seconds for the display to return to the last time you opened the bottle (c).

## DISCLAIMER

This e-pill Medication Reminder is only a supplement to any memory system that you currently use to remind yourself to take medications. The user maintains full responsibility to administer medication at the appropriate times. e-pill, LLC assumes no responsibility for any incidental or consequential damages, including but without limitation, damages or injuries resulting from malfunction of the product. e-pill, LLC's liability shall in no circumstances exceed the amounts paid by the user for the timer. Cap is not child-proof.