

# Instruction Manual

## e-pill ALRT Weekly Medication Reminder

### Easy Set-up

1. On the day and at the time of day you want to be reminded to take your weekly medication, remove pull tab from the back of the e-pill ALRT Weekly Medication Reminder.
2. The Reminder makes a 2-tone sound and then a single beep.
3. "7 days" will appear on the LCD and will commence countdown to 0 hrs. 7 days (168 hours) later your Weekly Medication Reminder will sound the alarm.



Size 6" x 5 1/2" x 1 3/4"

Your e-pill Weekly Medication Reminder is now set to go off ("alarm") 7 days from the time selected.

Your e-pill Medication Reminder can also be used as a Desk Clock and show current time (display will switch between number of days left and current time). In order for this feature to work, you must set the current time (see below).

### SET / ADJUST Desk Clock

- Press the "SET CLOCK" and the time displayed on the LCD will flash.
- Adjust the hour and the minute by pressing on the "HR" and the "MIN" buttons.
- After adjusting to the correct current time, press the "Set Clock" button again and the time is set. The LCD will stop to flash.

If the "Set Clock" button is not pressed a second time within 60 seconds after you have pressed the "HR" or "MIN" button, the flashing time displayed on the clock will be set as the current time.

### ACKNOWLEDGE ALARM (Silencing the Weekly Alarm)

- When the alarm sounds, take your weekly medication or treatment.
- Press "ALRT" button, this will silence the alarm.
- LCD will ONCE AGAIN DISPLAY 7 days or the number of days until the next alarm time (if you are late in acknowledging the alarm).

### CHANGE ALARM (Day and Time of the Weekly Alarm Setting)

To set the e-pill Weekly Reminder Alert to the right day and time of day, follow the this procedure (make sure the pull tab has been removed from the back of the e-pill device).

- At time and day desired to take medication every week – move the switch to "SET" position.
- Press and hold "ALRT" button until you hear a two-tone sound – then release.
- After releasing, a single beep will sound to confirm the Reminder is programmed.

- Move the switch back to the “LOCK” position. Your e-pill ALRT Weekly Medication Reminder has now been set to alarm 7 days (168 hours) from now.

The LCD display will show “7 Days” to the next alarm. If the current time has been set, the display will switch between number of days left and current time.

### **CLEAR ALARM (Erasing) the Weekly Alarm Setting**

- Move switch up to the “SET” position.
- Press and hold the “ALRT” button, continue to hold as you hear a single beep and then a two-tone sound.
- Now release the button.
- The Reminder is now cleared. Move the switch back to the “LOCK” position.

If the current time was programmed on the e-pill Weekly ALRT Reminder, it will continue to be displayed.

### **Additional Features**

**Alarm Duration / Missed Alarm:** If you do not hear your Weekly Reminder when it sounds, the LCD will continue to blink “0 hrs” indefinitely and the alarm will sound the alert for 30 seconds every 10 minutes (every 15 minutes after the first hour).

This will continue until you press the “ALRT” button, which will sound the alarm again. Then press the “ALRT” button again to acknowledge the event (see above).

**Battery:** Batteries are included and are replaceable. Battery life will vary depending on usage but with a fresh battery will exceed 1 year. Your e-pill ALRT Weekly Medication Reminder contains one Lithium Battery CR2032, or equivalent. To replace the battery, loosen the small center screw on the back of Reminder and remove the small cover. Remove the old battery and insert a new CR2032, with the negative (-) side facing down and the positive (+) side facing out; then replace the cover and tighten the screw. The Reminder will beep and the 7 Days will show in the LCD and the countdown begins. If this is not the desirable time for your alert, follow the instructions above for CHANGE ALARM.

**Warranty:** This e-pill ALRT Weekly Medication Reminder (excluding batteries) is backed by a 12 MONTH WARRANTY by the manufacturer. In the unlikely event that your ALRT Weekly malfunctions, during the first 12 months of use, we will replace your e-pill ALRT Weekly Medication Reminder with a new one.

*The e-pill ALRT Medication Reminder should be used as a supplement for remembering to take or administer medication on time, and should not be relied upon solely.*