

BlueMax
HD Lighting™

70w Desk Lamp

Another Great Product Manufactured By
Full Spectrum Solutions, Inc.



- Produces 10,000 Lux for light therapy.
- Fully dimmable—Light equivalent to 300w incandescent.
- High Definition lighting for visual acuity.
- Energy-saving bulb lasts up to 10,000 hours.

Full Spectrum Solutions, Inc.
www.fullspectrum solutions.com

Returns - Returns for exchange or refund must be handled by the store where the lamp was purchased.

Limited Lifetime Warranty:

This product is guaranteed to be free from defects in material and workmanship for life. Warranty on the bulb is one (1) year.

Full Spectrum Solutions, Inc. will replace any defective parts within the lifetime of the original purchaser under normal indoor use. This warranty covers all original labor, should you send the lamp into us for repairs, and parts excluding the light bulbs. The warranty does not cover damage to the unit caused by abuse, faulty household wiring, weather related events, or conditions related to normal wear. All warranty parts must be shipped pre-paid **in the original packaging** and insured by the consumer. Any damage incurred during return shipping is the responsibility of the customer and will void any warranty claim. All shipping costs will be the responsibility of the purchaser for shipping to and from the manufacturer. All covered labor must be authorized in advance by the Full Spectrum Solutions, Inc. Return and Warranty Department.

If your bulb arrives broken, you have two weeks from the shipment date to make a damage claim and you must contact the company that shipped the lamp to you. Bulbs broken by the customer are not covered under warranty.

This warranty only covers failures as stated above, which occur during normal use and in accordance with application guidelines. Full Spectrum Solutions, Inc. shall not be held liable for incidental or consequential damages resulting from the use of this product or arising out of any breach of this warranty. Full Spectrum Solutions, Inc. makes no other representation or warranties, expressed or implied and makes no warranty of merchantability or fitness for any particular purpose. We make no medical claims regarding use of this unit. This warranty is non-transferable. Purchaser agrees by the act of purchasing this product to the terms stated herein.

Dear Valued BlueMax™ Customer,

Congratulations! You have just purchased a quality BlueMax™ High Definition full spectrum lamp from Full Spectrum Solutions, Inc. Your lamp has been produced with the finest components available.

Be sure to retain your receipt and original packaging for warranty purposes. Put your receipt with this instruction manual and keep it accessible should you need it for future reference. If you have a warranty claim or you have questions about the product, please contact the dealer where you purchased your light.

INCLUDED PARTS:

- Desk lamp fixture
- Light bulb

ASSEMBLY AND START UP:

- Allow the lamp and bulb to reach room temperature before assembly.
- Do not plug the lamp into an electrical socket until fully assembled.
- Remove all the plastic from the lamp and unwrap the bulb.
- When installing the bulb, hold the bulb close to the plastic base rather than the glass end and push it into the head of the lamp. Make sure that the bulb is in all the way and the end of the bulb is secured by the metal clip. **See diagram on page 2.**
- Your lamp is now assembled.
- Insert the plug into desired electrical socket and turn the switch to the highest setting (largest dot).
- Leave the lamp on the highest setting for two hours to “burn in” the bulb properly for dimming.
***Note** - You only need to do this when you install a new bulb.
- After you have finished the two-hour “burn in” period you can start dimming the lamp.
***Note** – The white silicone substance between the bulb tubes acts as a stabilizer and should not be removed.

General Recommended Use for Light Therapy:

- Light therapy needs vary from person to person. Generally, 20-40 minutes of therapy in the morning is recommended. You do not need to stare directly at the lamp. To experience light therapy, you will want it in your peripheral while reading, enjoying a meal, etc.
- BlueMax™ 70w lamps produce 10,000 Lux at a distance of 18” and 5,000 Lux at 25”. If you find it uncomfortable to use the light on high, simply turn the light down or sit further away.

***Note** – This may require you to increase your light therapy session.

- You may use your lamp during any time of the day. However, it is recommended that you reduce the light to lower levels in the evening if you have been exposed to higher levels of light during the day.
- Using your light for work projects, displays and art, as well as light therapy should be handled like any other lighting environment application.

*** Important Notice** – We make no claims as to the medical benefits of these units and recommend that you contact a medical professional for specific recommendations and use. Our recommendations are based on user averages and can vary between users. It is your responsibility to determine what is best for you.

Bulb Installation and Replacement:

Make sure to unplug the lamp before installing or replacing the bulb.

Place the base of the bulb into the socket carefully. Adjust the metal clamp around the bulb. Use your fingers to pull the clamp around the bulb. To complete installation and ensure a firm connection press down on the plastic base of the bulb until it snaps in. You should feel or hear it snap into place. **Refer to diagram below.**

To Replace:

To change the bulb, carefully release the bulb from the metal clips while holding the bulb at the plastic base. Slide the bulb away from the lamp socket.



Troubleshooting and Reset Instructions:

Step 1: Check Bulb for Damage

Make sure that the bulb does not have any visible damage, such as a crack dark spots or a hole.

Step 2: Check Bulb Installation

Take a close look at the bulb and the lamp holder to check the connection. The male connector at the base of the bulb must be inserted so that you cannot see any part of it. You should feel or hear it snap into place. **Refer to diagram on page 2.**

If you cannot get the bulb to stay in the socket, see if the metal clip that holds the bulb is off-center. It can be repositioned if it is. The metal clip is screwed onto a metal strip that runs the length of the head of the lamp. You can remove the metal strip from the back of the head of the lamp by pulling down on the metal clip. The material that is on the back of that metal strip is a silicone substance so it should re-adhere after you move the metal strip to re-center it. The silicone substance is only in place to prevent the metal strip from rubbing on the head of the lamp during transit.

Step 3: Reset Lamp

Make sure the dimmer switch is in the “off” position. Unplug the lamp and let it sit for five (5) minutes. Turn the switch to the highest setting (largest dot), then plug in your lamp. The light should come on immediately.

Still Having an Issue?

Full Spectrum Solutions can be reached at (888) 574-7014 or by email at warranty@fullspectrumolutions.com. In order to better assist you, please be prepared to provide the model number, and the name of the store where the lamp was purchased.

Warranty Claims, Damage Claims, and Returns

If you purchased this product from a dealer other than Full Spectrum Solutions, Inc., you must contact that dealer.

If you purchased from Full Spectrum Solutions, Inc.:

Go to www.fullspectrumolutions.com for details on how warranty claims, damage claims, and returns are handled, or email: warranty@fullspectrumolutions.com.

Warranty Claims - IF the troubleshooting and reset instructions on page 3 failed to resolve your issue, contact the manufacturer, Full Spectrum Solutions, Inc. For warranty service, email warranty@fullspectrumolutions.com.

Damage Claims - Damage claims must be handled by the store where the lamp was purchased. We must receive all warranty parts shipped pre-paid by the customer. A Return Authorization Number must be obtained from Full Spectrum Solutions prior to returning product.