


# PREPSTER™

## Electrode Skin Prep Pad



Use **PREPSTER** to gently abrade patient skin prior to ECG, EEG, Audiology, Stress, TENS, Sleep and other studies

**PREPSTER** is a specially formulated skin prep pad designed for healthcare providers, Made in U.S.A. from Hypoallergenic material - **Latex-free**

Pre-test preparation with **PREPSTER** produces no mess, leaves no gritty residue and needs no clean-up

*“ Better Skin Preparation Means A More Accurate Reading ”*

“Preparing the skin to create Clean Contact has been shown to dramatically improve electrode test readings.”

“Exfoliation removes the unwanted, nonconductive skin cells, with minimal time and effort... several gentle strokes can effectively decrease impedance...”

“Poor skin preparation leads to unwanted artifact and needless retesting. Proper skin prep techniques can improve patient assessment”

### **PREPSTER™** Skin Prep Pads

Convenient 2.5 inch x 1.75 inch size

Disposable single use pads, non-sterile

Prepster side backed with soft material, white color

80 **PREPSTERS** per zip-lock bag

**Made in U.S.A.**



Distributed by:

**medigear.com**  
Natick, MA USA